

# Statistical Data Analysis Project

Alexandre De Cuyper

University of A Coruña

February 22, 2024

# Introduction

- ▶ Analyze the powder X-ray diffraction data from Complex Materials,
- ▶ Objectives and hypotheses: Investigating if there are significant differences in peak positions and intensity among crystals,
- ▶ Importance of the analysis: Understanding the change in structure of the crystals when changing the composition.

# Data Processing

- ▶ Exporting the data from MATCH!
- ▶ Importing the data with R.


```
1 > head(data)
2   2_theta Ni75Co25 Ni50Co50 Ni25Co75   Co
3 1    8.00    20.0    18.00    24.00 19.0
4 2    8.05    19.5    21.98    23.01 20.5
5 3    8.10    19.0    26.00    22.00 22.0
6 4    8.15    19.0    27.49    21.50 25.0
7 5    8.20    19.0    29.00    21.00 28.0
8 6    8.25    17.0    24.53    24.98 24.5
9
```

# Exploratory Data Analysis

- ▶ Visualizations and summary statistics.
- ▶ Identification of patterns and trends.


```
1  # Exploratory Data Analysis (EDA)
2  ggplot(data, aes(x = '2_theta')) +
3  geom_ribbon(aes(ymin = 0, ymax = 'Ni25Co75', fill = "Ni25Co75"),
4  alpha = 0.1, color = "green"
5  ) +
6  geom_ribbon(aes(ymin = 0, ymax = 'Ni50Co50', fill = "Ni50Co50"),
7  alpha = 0.1, color = "red"
8  ) +
9  geom_ribbon(aes(ymin = 0, ymax = 'Ni75Co25', fill = "Ni75Co25"),
10 alpha = 0.1, color = "blue"
11 ) +
12 geom_ribbon(aes(ymin = 0, ymax = 'Co', fill = "Co"),
13 alpha = 0.1, color = "yellow"
14 ) +
15 labs(
16   title = "Powder X-ray Diffraction of Ni_(1-x)Co_x perovskites",
17   x = "2_theta",
18   y = "Intensity",
19   fill = "Composition"
20 ) +
21 scale_fill_manual(values = c(
22   "Ni75Co25" = "blue", "Ni50Co50" = "red",
23   "Ni25Co75" = "green", "Co" = "yellow"
24 )) +
25 theme_minimal()
26
```

# Visualization of the Data



# Visualization of the Data

Powder X-Ray Diffraction of  $\text{Ni}_{(1-x)}\text{Co}_x$  perovskites




# Exploratory Data Analysis

- ▶ Visualizations and summary statistics.
- ▶ Identification of patterns and trends.

```
1 #Threshold value
2 +
3 geom_hline(
4   yintercept = threshold,
5   linetype = "dashed", color = "black"
6 )
7
```

# Visualization of the Data






# Peak Identification

```
1 find_peaks <- function(data) {
2   if (length(data) < 2) {
3     return(NULL) # No peak in lists with 0 or 1 element
4   }
5   peaks <- rep(0, length(data))
6   # Create a vector to store the value and index of the peaks
7   for (i in 2:(length(data) - 1)) {
8     if (data[i] > data[i - 1] && data[i] > data[i + 1]) { # Looking for a peak
9       peaks[i] <- data[i]
10    }
11  }
12  # Checking if the first and last value is a peak or not
13  if (data[1] > data[2]) {
14    peaks[1] <- data[1]
15  }
16  if (tail(data, 1) > tail(data, 2)[1]) {
17    peaks[length(peaks)] <- data[length(data)]
18  }
19  return(peaks)
20 }
```

# Peak Identification Results




# Clustering Process

- K-means clusters data by iteratively assigning points to the nearest cluster center and updating those centers until convergence, aiming to group similar points together.

```
1
2 # Specify the number of clusters (k)
3     k <- 22
4
5 # Perform k-means clustering based only on 2_theta
6     cluster_assignments <- kmeans(data_for_clustering,
7         centers = k,
8         nstart = 4
9     )$cluster
10
11
```


# Clustering Results

- Overview of the clustering results.
- Interpretation of clusters and their characteristics.




# Clustering Results

- Overview of the clustering results.
- Interpretation of clusters and their characteristics.



# Clustering Results

- Overview of the clustering results.
- Interpretation of clusters and their characteristics.



# ANOVA Test

```
1 # Perform ANOVA
2 anova_result_composition <- aov('2_theta' ~ Composition *
   cluster, data = non_zero_data_long)
3
4 # Print ANOVA summary
5 summary(anova_result_composition)
6
```

## ANOVA Results

Source	Df	Sum Sq	Mean Sq	F Value	Pr(>F)
Composition	3	156.5	52.15	81.746	0.000114***
Cluster	15	2896.1	193.08	302.627	2.33e-06***
Composition:cluster	31	0.9	0.03	0.044	1.000000
Residuals	5	3.2	0.64		


**Table:** ANOVA Results for 2theta with Composition and Cluster.



## Analysis of ANOVA Results


- ▶ The Composition factor shows a significant effect on the response variable (p-value = 0.000114).
- ▶ The Cluster factor has a highly significant impact on the response variable (p-value = 2.33e-06).
- ▶ The interaction between Composition and Cluster is not significant (p-value = 1.0000).
- ▶ Residuals indicate the unexplained variance in the model.

# Results



# Peak Intensity

- ▶ Powder X-ray Diffraction (PXRD) data often contains peaks that correspond to specific crystallographic planes.
- ▶ Peak intensity is a crucial parameter in PXRD analysis, reflecting the abundance or concentration of particular crystallographic phases.



## Peak Intensity Analysis

Composition	Mean Intensity	Variance Intensity
Co	99.8	3653.0
Ni25Co75	116.0	8636.0
Ni50Co50	104.0	8514.0
Ni75Co25	109.0	6062.0

Table: Mean and Variance Intensities for Each Composition.

# ANOVA Results

- ▶ Analysis of Variance (ANOVA) was performed to assess the impact of 'Composition' and 'cluster' on the 'Intensity' variable.
- ▶ Statistical significance was evaluated based on p-values.

```
1 > anova_result <- aov(Intensity ~ Composition * cluster, data = non_zero_data_long)
2
3
4 > summary(anova_result)
```

Factor	Df	Sum Sq	Mean Sq	F value	Pr(>F)
Composition	3	9496	3165	14.80	0.012440*
Cluster	9	220856	24540	114.78	0.000182***
Composition:Cluster	21	31651	1507	7.05	0.035380*
Residuals	4	855	214	—	—

Table: ANOVA Summary

## ANOVA Results

- ▶ **Composition:** The p-value (0.012440) indicates a significant difference in means across 'Composition' levels.
- ▶ **Cluster:** A very low p-value (0.000182) suggests significant differences in means across clusters.
- ▶ **Interaction:** The interaction between 'Composition' and 'cluster' is significant (p-value = 0.035380).

# Conclusion

## Summary of key findings.

- ▶ Significant peak position and intensity variation depending on the composition,
- ▶ Implies a change in the crystalline structure of the material.

## Limitations and areas for future research.

- ▶ **Ni** composition for further analysis,
- ▶ Improve the peak finding algorithm,
- ▶ Random variability in the clustering process.

Any Questions?